

Trends in Cochlear Implant Complications and Clinical Outcomes from 2016 to 2024

Daniel R. S. Habib, BA; Anthony E. Bishay, BS; Ankita Patro, MD;
Alexander J. Langerman, MD SM FACS; Aaron C. Moberly, MD;
David S. Haynes, MD MMHC FACS; Kareem O. Tawfik, MD

Daniel Roy Sadek Habib

MS4 at Vanderbilt University School of Medicine

✉ Daniel.r.habib@Vanderbilt.edu

✕ @danielrshabib



DISCLOSURE INFORMATION

- Standards for Integrity and Independence in Accredited Continuing Education



The AAO-HNSF is committed to creating high-quality education that is independent of industry influence. All persons who have been in control of this educational content for this accredited CE activity has been asked to disclose all financial relationships with any ineligible companies they have had over the past 24 months.

The ACCME defines ineligible companies as those whose primary business is producing, marketing, selling, re-selling, or distributing healthcare products used by or on patients. There is no minimum threshold; all financial relationships, regardless of amount, with ineligible companies.

The ACCME also requires that AAO-HNSF manage any reported conflict and eliminate the potential for bias during the educational activity. All of the relevant financial relationships listed for these individuals have been mitigated. The disclosure information is intended to identify any commercial relationships and allow learners to form their own judgments. However, if you perceive a bias during the educational activity, please report it on the evaluation.

***Kareem O. Tawfik, MD**, faculty for this accredited education activity has no relevant financial relationships with ineligible companies to disclose.*



Cochlear Implantation Prevalence and Complications



- 180,000+ total cochlear implants (CIs) in U.S.^{1,2}



- Single-center CI complication rates: 9-20%³⁻⁶



Cochlear Implant Adverse Events in MAUDE



- Before 1998 vs in 2002: ↓ device failures, ↑ infection⁸
- 2000–2010: ↑ idiopathic adverse events⁹
- 2010–2020: Infection, extrusion, and neurologic complications differed by manufacturer.¹⁰



Voluntary Field Corrective Action of HiRes Ultra and Ultra 3D

When a company announces a recall, market withdrawal, or safety alert, the FDA posts the company's announcement as a public service. FDA does not endorse either the product or the company.

[Read Announcement](#)

Summary

Company Announcement Date:	February 20, 2020
FDA Publish Date:	March 10, 2020
Product Type:	Medical Devices
Reason for Announcement:	Reports of hearing performance degradation
Company Name:	Advanced Bionics
Brand Name:	HiRes Ultra and Ultra 3D
Product Description:	Cochlear implants

Research Question

How did U.S. CI complications and clinical outcomes differ before and after a major VFCA in 2020?



Hypothesis

After a major VFCA in 2020, there will be decreases in U.S. CI failures and required interventions.



Methods



Manufacturer and User Facility Device Experience
(MAUDE) Database

Data

- Filter by CI product codes
[MCM & PGQ] from 2016–2024

Methods



Manufacturer and User Facility Device Experience
(MAUDE) Database

Data

- Filter by CI product codes [MCM & PGQ] from 2016–2024

Analysis

- Quasi-Poisson regressions for patterns in raw counts [2016–2020 & 2020–2024]

Methods



Manufacturer and User Facility Device Experience
(MAUDE) Database

Data

- Filter by CI product codes [MCM & PGQ] from 2016–2024


Analysis

- Quasi-Poisson regressions for patterns in raw counts [2016–2020 & 2020-2024]

Outcomes

- 1) Patient Complications
- 2) Device Complications
- 3) Post-Op Clinical Outcomes
 - Required Intervention
 - Hospitalization
 - Life-Threatening Event/Disability/Death
 - Other

Number of Top Device Complications




- 33,487 device complications
- 2016–20: ↑ Fit/malposition/detachment, Electrical, Handling/software issues
- 2020–24: ↑ Fit/malposition/detachment; ↓ Unexpected result, Quality/safety issues

Device Complication	2016-20	2020-24
	IRR (95% CI), P-Value	IRR (95% CI), P-Value
Unexpected Result/Output	1.03 (0.92-1.14), p=.645	0.69 (0.58-0.82), p=.024
Fit/Malposition/Detachment	1.55 (1.39-1.74), p=.005	1.08 (1.06-1.11), p<.001
Electrical	1.61 (1.34-1.93), p=.014	1.25 (1.08-1.44), p=.057
Handling/Software	2.34 (1.91-2.87), p=.004	0.99 (0.79-1.24), p=.918
Quality/Safety/Nonstandard Use	1.06 (0.90-1.24), p=.534	0.83 (0.74-0.93), p=.047



Number of Top Patient Complications




- 34,426 patient complications
- 2016–20: ↑ Device failure/unexpected result
- 2020–24: ↑ Healing/skin issues

Patient Complication	2016-20	2020-24
	IRR (95% CI), P-Value	IRR (95% CI), P-Value
Device Failure/Unexpected Result	1.15 (1.06-1.24), p=.042	0.92 (0.88-0.97), p=.058
Healing/Skin	1.27 (1.07-1.52), p=.075	1.22 (1.08-1.38), p=.047



Number of Clinical Outcomes



- 26,377 clinical outcomes
- 2016–20: ↑ Required intervention, Other
- 2020–24: ↑ Hospitalizations, Other

Outcome	2016-20	2020-24
	IRR (95% CI), P-Value	IRR (95% CI), P-Value
Required Intervention	1.10 (1.06-1.15), p=.022	1.00 (0.97-1.02), p=.747
Other	1.12 (1.03-1.21), p=.006	1.32 (1.14-1.55), p=.038
Hospitalization	1.62 (0.55-4.73), p=.445	1.50 (1.14-1.99), p=.064



Limitations



- Retrospective
- Reporting bias (96.6% mandatory¹¹)
- Missing surgical/comorbidity data



Conclusion



- Stable total number of complications and required interventions



- ↑ Fit/malposition/detachment, Healing/skin issues, Hospitalizations



- ↓ Unexpected device result/output



References

1. Nassiri AM, Marinelli JP, Lohse CM, Carlson ML. Incidence of Cochlear Implantation Among Adult Candidates in the United States. *Otol Neurotol*. 2023;44(6):549-554. doi:10.1097/MAO.0000000000003894
2. National Institute on Deafness and Other Communication Disorders. Cochlear Implants. Published 2024. Accessed July 27, 2025. <https://www.nidcd.nih.gov/health/cochlear-implants>
3. Achiques MT, Morant A, Muñoz N, et al. Complicaciones y fallos de la implantación coclear. *Acta Otorrinolaringológica Española*. 2010;61(6):412-417. doi:10.1016/j.otorri.2010.07.005
4. Farinetti A, Ben Gharbia D, Mancini J, Roman S, Nicollas R, Triglia J-M. Cochlear implant complications in 403 patients: Comparative study of adults and children and review of the literature. *Eur Ann Otorhinolaryngol Head Neck Dis*. 2014;131(3):177-182. doi:10.1016/j.anorl.2013.05.005
5. Chen DS, Clarrett DM, Li L, Bowditch SP, Niparko JK, Lin FR. Cochlear Implantation in Older Adults. *Otol Neurotol*. 2013;34(7):1272-1277. doi:10.1097/MAO.0b013e3182936bb2
6. Halawani R, Aldhafeeri A, Alajlan S, Alzhrani F. Complications of post-cochlear implantation in 1027 adults and children. *Ann Saudi Med*. 2019;39(2):77-81. doi:10.5144/0256-4947.2019.77
7. U.S. Food & Drug Administration. Manufacturer and User Facility Device Experience (MAUDE) Database. Published 2024. Accessed February 2, 2025. <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfmaude/search.cfm>
8. Tambyraja RR, Gutman MA, Megerian CA. Cochlear Implant Complications. *Arch Otolaryngol Neck Surg*. 2005;131(3):245. doi:10.1001/archotol.131.3.245
9. Causon A, Verschuur C, Newman TA. Trends in Cochlear Implant Complications. *Otol Neurotol*. 2013;34(2):259-265. doi:10.1097/MAO.0b013e31827d0943
10. Jinka S, Wase S, Jeyakumar A. Complications of cochlear implants: a MAUDE database study. *J Laryngol Otol*. 2023;137(11):1267-1271. doi:10.1017/S0022215123000828
11. Kavanagh KT, Brown Jr RE, Kraman SS, Calderon LE, Kavanagh SP. Reporter's occupation and source of adverse device event reports contained in the FDA's MAUDE database. *Patient Relat Outcome Meas*. 2019;Volume 10:205-208. doi:10.2147/PROM.S212991



Thank you!

Questions?

Daniel Roy Sadek Habib

MS4 at Vanderbilt University School of Medicine

✉ Daniel.r.habib@Vanderbilt.edu

✕ @danielrshabib

Slides



SCAN ME

Advanced Bionics Complications

